

**ARTG LISTED**  
**372335**

# COVID-19 Antigen Nasal Test Kit (Nasal)

**Non-Intrusive**



## Benefits

- User-friendly, no equipment required
- Fast results within 15 minutes

## Clinical Evaluation

- Relative sensitivity: 98.1% (93.4%~ 99.5%)
- Relative Specificity: 99.8% (98.6%~100.00%)
- Overall Agreement: 99.4% (98.3%~99.8%)  
\*95% Confidence Interval

**Table 1 Clinical Evaluation**

	RT-PCR		Total
	Positive	Negative	
Positive	104	1	105
Negative	2	401	403
Total	106	402	508

## Test Procedure

