



P2 RESPIRATOR WITH EARLOOPS



TECHNICAL SPECIFICATION SHEET

Key Features

- Filters >98% of airborne particles
- Level 2 fluid resistance
- Tested and certified in Australian laboratories
- Comfortable earloops and adjustable nose clip for a snug fit
- Individually wrapped for sharing and on-the-go


Material Composition


- Layers: Polypropylene (PP) and cotton
- Earloops: Spandex
- Nose clip: PP coated twin wire
- Free from latex


Country of Origin

- Made in Australia from imported materials

Directions for Use

- Step 1
 After washing hands, remove mask from packaging.

- Step 2
 Press mask against face and place earloops over ears.

- Step 3
 Push nose guard securely down for a tight fit.



Approvals and Standards

- Certified P2 Particulate Respirator
- Meets: AS/NZS 1716:2012 (Equiv. EN149:2001, GB2626-2019)
- Meets: ASTM F1862/F1862M-12/ISO22609
- Meets: EN14683:2014, Appendix C
- TGA number: ARTG 352290
- CertMark number: TT10061





P2 RESPIRATOR WITH EARLOOPS



Packaging Information

- 25 Masks (Retail Box)
- 25 masks per box
- 600 masks per shipper
- 12,000 masks per pallet
- 200,400 masks per 20' container
- 450,000 masks per 40' container

- Mask weight: 4.6g
- Box weight: 200g
- Shipper weight: 4.8kg
- Pallet weight: 96 kg

100 Masks (Wholesale Box)

- 100 masks per box
- 285,600 masks per 20' container
- 612,000 masks per 40' container

Shelf Life and Storage

- Store in a cool dry place, away from direct sunlight.
- Minimum 2 years to maximum 5 years shelf life from date of manufacture
- Refer date of manufacture (DOM) on base of packaging



6 Precision Lane, Notting Hill, VIC 3168
<https://www.ausmedhealth.com.au/>
contact@ausmedhealth.com.au

