

Maxi Clear¹⁵

15 DRUG URINE TEST

The future in Australian drug testing

The Maxi Clear 15 will detect 15 drugs of abuse plus 6 adulterants. Certified to Australian Standards AS/NZS4308:2008 requirements.

Drugs tested

Amphetamine	Methamphetamine
Cocaine	THC Marijuana
Opiates	Ecstasy
Benzodiazepines	Fentanyl
Ketamine	Oxycodone
Methadone	Buprenorphine
Barbiturates	Phencyclidine
Propoxyphene	



Maxi Clear 15 Urine Drug Test Kit is a comprehensive drug testing kit for all applications. A reliable onsite test with results within minutes. The larger collection cup makes sampling more convenient for the donor, and the cup includes specimen validity tests to detect specimen adulteration.

Benefits

- Simple, easy to use, accurate urine drug test cup
- Fast one-step operation with results in minutes
- Larger cup for easier sample collection
- Built-in specimen validity test to detect Specimen Adulteration Oxidants, Specific Gravity, pH, Nitrite, Glutaraldehyde, Creatinine
- Excellent pricing and quality

Ideal for

- Pre-employment drug screening
- Workplace drug testing programs
- Random drug tests
- Medical clinic drug testing
- Reasonable suspicion drug testing
- Substance abuse screening in hospitals

With the additional synthetic drugs now becoming more prevalent in Australia, the Maxi Clear 15 Drug Urine test cup is a complete drug testing kit. A simple one step immunoassay for detection of drugs and drug metabolites in urine.



Interested in a kit with a different drug combination?

Our range of tests are developed to meet current market demands. If you are interested in discussing a customised kit specific to your drug testing needs, please let us know.