




Determine your Pacey Cuff Size

STEP 1



Circumference

Measurement is done mid-shaft around a non-erect (flaccid) penis.


Penis Circumference

SIZE	INCHES	CENTIMETERS
S	1 - 2.75"	3 - 7 cm
M	2.75 - 4	7 - 10 cm
L	4 - 5.5"	10 - 14 cm

* if your measurement is between two sizes we recommend choosing the larger size.

Ensure you have the length required


STEP 2



Measurement is from base to tip of non-erect penis.

Minimum Penis Length

SIZE	INCHES	CENTIMETERS
S	1.25"	3 cm
M	1.75"	4 cm
L	2"	4.5 cm



Use this strip to accurately measure your circumference