



## SMART WEARABLE THERMOMETER

Fever monitoring wherever you are

### Remote & continuous monitoring

In the form of a soft, comfortable and supple patch, Tucky is placed under the armpit and allows to monitor evolution of body temperature.

Real time temperature curve is displayed on a remote smartphone app.

An alert is raised when fever exceeds a threshold temperature.



### No distance limit

Fever can be monitored either at home (5-10m though Bluetooth) or remotely (at work) via internet on a 2nd smartphone (secured cloud data synchronization).



### An intuitive app

- Multi profiles
- Multi devices
- Data history
- Curve / timeline views
- Health book...

## A secured medical device



### Class IIA Medical Device

CE marking / ISO13485 Certified



### Patented Technology




Ultra low wave emission



### Secured data

Certified Healthcare Host

## The most complete thermometer of the market

Thermomètres	<b>+TUCKY</b>	Ear / Temporal	Traditional
			
24h continuous monitoring	✓	✗	✗
Remote tracking	✓	✗	✗
Fever alert	✓	✓	✗
Temperature Color code	✓	✓	✗
Data history	✓	✓	✗
Position monitoring	✓	✗	✗

- Flexible, comfortable, hypoallergenic
- Medical certification, accuracy 0,1°C
- Ultra low wave emission (1000 lower than a smartphone)
- Rechargeable (autonomy > 5 day)

## + product

Tucky is also a **position monitor**. It warns when newborn is no longer in the recommended sleeping position.

## Tucky box & adhesive refill kit

### Tucky box includes:

- 1 thermometer
- 5 adhesives
- 1 USB recharge cable
- 1 manual

**Adhesive refill kit:** x15 (ref=ADH21)



Model: TUCKY-21W – Class IIA medical device  
 Manufactured by: e-TakesCare SAS, 35 rue des Chantiers, 78000 Versailles, France  
 More information: [www.e-takescare.com](http://www.e-takescare.com)

