

RESTLEG STRAP FOR WHEELCHAIR



DESCRIPTION

The restleg strap for wheelchair / chair offers rest for the area of the calves.

MODEL / REFERENCE

307201

COLOURS

Blue

EFFECTS AND ADVANTAGES

It offers a firm and comfortable support in the area of the calves with security lock on the back.

INDICATIONS OF USE

To place the article correctly, the patient must be placed comfortably in the wheelchair / chair and the legrest band must pass behind the calves adjusted to the most comfortable height for the patient.

COMPOSITION

Padded fabric: 20% Antelina (100% polyester), 20% Velur (100% polyamide) and 60% foam.

Breathable: 100% polyester

Straps: 100% Polypropylene

MAINTENANCE

Simply wash it periodically by hand or machine, at a maximum temperature of 40°C. Dry horizontally and without heat sources.