









Surgical Face Mask Non-Sterile

Description: This Surgical Mask provides effective respiratory protection against solid and liquid aerosols that can contain oil. They are intended to protect from the transfer of microorganisms, body fluids and particulate material.

FEATURES:

-  3-layer construction
-  Two layers of polypropylene non-Woven fabric with a layer of polypropylene melt-blown filter cloth in the middle
-  The mask is held in place over the user's mouth and nose by two elastic ear loops
-  The adjustable wire and plastic nosepiece helps to form a good seal and comfortable wear
-  The disposable masks are designed to provide you with maximum comfort
-  Manufactured in an ISO 13485 certified facility

TECHNICAL SPECIFICATIONS:

| | |
|-----------------------|----------------------|
| ASTM Level | 2 |
| Outer Layer Material | SPP |
| Inner Layer Material | SPP |
| Basis Weight (gms) | 73 |
| Color | Blue |
| Middle Layer Material | Meltblown |
| Nose Wire | Wire with PE Plastic |
| Ear Loop | Spandex |