












Star 6" Adjustable Rollator

Features:

- Quick-Fit adjustable design
- 6 holes handle adjustment
- Easy to pack up for travel
- 3 seat height adjustment
- Locking pin to keep walker folded
- Aust. STD Tested 493360



Length (cm)		65.5-73.1	Carton Size (cm)	37*21*16	Front Castor Diameter (cm)		15.2/6"	
Width (cm)		61.7	Seat Depth (cm)		34.8	Rear Castor Diameter (cm)		15.2/6"
Max Height (cm)		93.5	Seat Width (cm)		37	Max Load (kg)		150
Minimum Height (cm)		81.5	Seat Height from the Ground (cm)		49.3-57	Net Weight (kg)		7.3

Star 8" Adjustable Rollator

Features:

Quick-Fit adjustable design

6 holes handle adjustment












Easy to pack up for travel

3 seat height adjustment

Locking pin to keep walker folded

Aust. STD Tested 493360



Length (cm)		66.8-71.9	Carton Size (cm)	37*16*21	Front Castor Diameter (cm)		20/8"	
Width (cm)		61.8	Seat Depth (cm)		34.8	Rear Castor Diameter (cm)		20/8"
Max Height (cm)		96.8	Seat Width (cm)		36	Max Load (kg)		150
Minimum Height (cm)		80.7	Seat Height from the Ground (cm)		50.8-55.8	Net Weight (kg)		8.3

Star 8" HD Adjustable Rollator

Features:

Quick-Fit adjustable design

6 holes handle adjustment












Easy to pack up for travel

3 seat height adjustment

Locking pin to keep walker folded

Aust. STD Tested 493360



Length (cm)		66.8-72	Carton Size (cm)	43*23*16	Front Castor Diameter (cm)		20/8"	
Width (cm)		70.8	Seat Depth (cm)		38	Rear Castor Diameter (cm)		20/8"
Max Height (cm)		96.8	Seat Width (cm)		46	Max Load (kg)		204
Minimum Height (cm)		80.8	Seat Height from the Ground (cm)		50.8-55.8	Net Weight (kg)		9.5