

## Walking Aids



### Forearm Crutch

- \* Lightweight aluminium
- \* Ergonomic handle
- \* Double adjustment, length & cuff
- \* Non-slip rubber ferrule
- \* Max weight user: 110kg

Standard: 86-116cm  
Tall: 108-136cm

### Axillary Crutches

- \* Lightweight aluminium
- \* Adjustable hand support
- \* Locking pin adjustment
- \* Non-slip rubber ferrule
- \* Max user weight: 100kg



#### Product Code:

**CRAXS** Small 94-114cm

**CRAXM** Med 114-134cm

**CRAXL** Lge 134-154cm



### Bariatric Elbow Crutches

- \* Lightweight aluminium
  - \* Wedge shape PVC handle
  - \* Closed cuff
  - \* 14 height adjustments
  - \* 3 cuff adjustments
  - \* Non-slip rubber ferrule
  - \* Max weight user: 222kg
- Handgrip height 635-915mm  
Overall Height 522-1210mm