

Bathroom

Rotating Seat Shower Chair

This Shower Chair features a rotating seat which rotates in both directions 360 degrees with lock positions at 90 degree intervals making it highly recommended for being user friendly to both the user and the carer. Constructed from a durable oval shaped aluminium frame with a detachable back and arm rests, height adjustable extension legs and suction rubber tips to resist slips.

- * Oval shaped aluminium frame
- * 360 degree rotating seat
- * Drainage holes
- * Removable backrest
- * Removable armrests
- * Height adjustable extension legs
- * Non-slip stoppers
- * Weight capacity 136kg



Product Code: RM410

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
85-97.5	55	58	33	51	50-62	34.5	6	136

Measurements in cm, Weight in kg

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com