

## 2in1 Electric Rollator/ Wheelchair



Load Capacity	120kg
Width Between Armrest	54 cm
Widest Width	65 cm
Seat Height to floor	50 cm
Handle Height	780 – 880 mm
Seat size	44 x 30 cm (H x W)
Ground to seat	50 cm
Motor	2 x 250W Brushless
Battery	24V 6AH Litium ion



## 2in1 Electric Rollator/ Wheelchair

Battery Location	Seat Bag
Weight without battery	16.6 kg
Weight with Battery	18kg
Charging time	2-3 hrs
Driving Distance	Up to 20km Depending on weight
Max Speed	6 km/hr

---

Front Joystick	High efficiency Intelligent 360° rotationa
Rear Joystick	Touch sensor drive
Front & Rear Tyre	8" / 8" PU Solid
Climbing Slope	15°
Unfolded Size Wheel to wheel	700 x W 700 H 900 mm
Folded Wheel to wheel	700 x W 290 x H 900 mm
Footrest	YES Swing away
Brake Type	Regenerative Electromagnetic

---



## 2in1 Electric Rollator/ Wheelchair

---

Large Shopper Bag

YES

---

Material

Aluminium Alloy

---

Colours

Gold/Yellow – Blue

---

Warranty

12 months

---