



Adjusta Wedge Body Support

TGA Listed	Yes
Australian Register of Therapeutic Goods Number	No.174456 Device: Pillow
Length (of backrest portion)	70cm
Width	60cm
Size of Semi Circle Support	60cm wide x 70cm high x 18cm deep
Size of Triangular Support	60cm wide x 40cm high x 27cm deep
Inner	High-Density Open Cell Urethane Foam
Outer Cover	100% Cotton quilted cover
Country of Manufacture	Australian Made
Function	2 piece construction
Angle Option	15 degrees to 75 degrees angle
Outer Cover	Tailored Fitted Easy Wash Cotton Cover