



Back Huggar Chair Cushion - Lumbar & Lower Back Support Seat Cushion

Dimensions	Length 450mm X Height 338 mm X Width/depth 70 mm at thickest point Support column (extra) width 10 m
Inner	100% open cell 'breathing' urethane foam OR visco-elastic Memory Foam.
Dura-Fab Outer Cover	Woven polyester. It is a super-soft easy-clean washable fabric.
Steri-Plus Outer Cover	Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester backing.