

Bariatric Coccyx Cushion



Outer Cover - Dura-Fab

Woven polyester - Suede Feel

Outer Cover

Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester backing

Inner

Inner Foam Construction

2.5cm memory foam + 5cm high density max support traditional foam + 2.5cm memory foam

Large Size

50cm (w) x 45cm (d) x 10cm (h)

Extra Large Size

56cm (w) x 45cm (d) x 10cm (h)

Country of Manufacture

Australian Made