



Big & Strong Self-Propelled Wheelchair



Seat Width	55 cm
Seat Depth	42 cm
Seat Height	49 cm (without cushion)
Backrest Height (from ground)	43 cm (87 cm)
Overall Height	94 cm
Overall Width	74 cm
Overall Length (with footrests)	106 cm
Armrest Height (from ground)	68 cm
Folded Height	82 cm
Folded Width	31 cm
Folded Length	102 cm
Total Weight (with accessories)	16 kg
Rear Wheels	24 inch / 60 cm Puncture proof solid tire
Front Castors	19 cm Puncture proof solid tire
Material	Aluminium
SWL	175 kg
Included	Seat cushion Seatbelt Calf rest