



Breezy BasiX 2 Folding Wheelchair



wheelchair-type

Self Propelled

Seat Width

14", 16", 18", 20" (52cm)

Seat Depth

410 mm – 460 mm

Seat Height

425 mm – 510 mm

Seat Angle

0 degrees to 6 degrees

Back Rest Height

410 – 460 mm

Material

Aluminium

Seat Tilt

No

Max SWL

125 kg

brand

Sunrise Medical

manufacturers-part-number

074100-017

Frame Weight

8.4kg
