

C-101 - 3 Position Lift Chair



Weight/Height:	MAX Weight 147 kg (325 lbs)/163-183cm (5'4-6'0)		
Features:	Button back design		
Back Style:	Blown Fibre Button-Back		
Head/Arm Covers:	No	Back Width:	72 cm (28.5")
Seat Depth:	51 cm (20")	Seat Width:	51 cm (20")
Top of Back to Seat:	70 cm (27.5")	Seat-to-Floor:	47 cm (18.5")



C-101 Electric Recliner Lift Chair

Max SWL	147 kg
Seat Height	470 mm
Seat Depth	510 mm
Seat Width	510 mm
Overall Width	780 mm