

## Colibri Mobility Scooter



**Speed** 8 km/h

---

**Max Range** 16km

---

**Max climb angle** 6 degrees/10.5%

---

**Curb Climbing** 4.5cm

---

**Batteries** 2 x 18Ah (Batteries will install with the scooter before dispatch)

---

**Max SWL** 136kg

---

**Weight** 49kg (incl. batteries)

---

**Seat Width** 46.5cm

---

**Seat Depth** 40cm

---

**Seat Height** 38.5cm

---

**Overall Width** 50.5cm

---

**Overall Length** 101cm

---

**Warranty** 24 months

---