



Compact Electric Recliner Lift Chair



Max SWL 140 kg

Back Height Small: 680 mm / Medium: 740 mm /
Large: 800 mm

Seat Height Small: 480 mm / Medium: 520 mm /
Large: 520 mm

Seat Depth Small: 440 mm / Medium: 510 mm /
Large: 510 mm

Seat Width Small: 500 mm / Medium: 500 mm /
Large: 600 mm

Overall Width Small: 750 mm / Medium: 750 mm /
Large: 820mm

Overall Height 1060 mm

Brand KCare

Manufacturer's part number AC59018

Sizes Single Motor Small, Single Motor
Medium, Single Motor Large, Twin

**Motor Small, Twin Motor Medium,
Twin Motor Large**

Warranty

**A warranty of 5 years on the frame, 2
years on castors and motor 1 year on
fabric and struts.**
