



Days Seat Walker/Wheelchair – Transit 2 in 1



Max SWL 160kg (352.74 lbs)

Seat height 56 cm (22 inches)

Seat width 35 cm (13.7 inches)

Seat depth 33 cm (12.9 inches)

Overall height 80 – 93 cm (31.4-36.6 inches)

Overall width 58.5 cm (23 inches)

Overall depth 74 cm (29.1 inches)

Weight 8.8 kg (19.4 lbs)

Wheel 7.5"

Color Red

Warranty 12 months
