

## Diffuser Cushion - Pressure Diffusing Memory Foam Chair Cushion



<b>Medium Size</b>	44 cm x 42 cm x 7.5 cm
<b>Large Size</b>	49 cm x 46 cm x 7.5 cm
<b>Inner</b>	100% open cell 'breathing' urethane foam and visco-elastic Memory Foam.
<b>Dura-Fab Outer Cover</b>	Woven polyester. It is a super-soft easy-clean washable fabric.
<b>Steri-Plus Outer Cover</b>	Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester backing.