



Dinkum NAVIGATOR + T-100 Trailer



BRAND	Dinkum NAVIGATOR
Material	6061 Aerospace Aluminum Alloy
Unfolded Dimension(LxWxH)	850x600x935mm
Folded Dimension(LxWxH)	680x600x360mm
Width of Armrests	480mm
Seat Size	450x420mm
N.W.(Without battery)	26kg
Gross Weight	27.8kg
Max Speed	7km/h
Load Capacity	160kg
Motor	DC Brushless Motor 2x250W
Battery x 2	Lithium Battery 24V/10AH 240WH
Battery Weight	1.8kg
Driving Range	40km
Tyres	Front: 7" PU Rear: 12" PU
Ground Clearance	59mm
Turning Radius	Own circle
Grade Ability	10°
