



Dinkum R120 Reclining Electric Wheelchair



MODEL	Dinkum R120Plus – Reclining Electric Wheelchair
Motor	2 x 250W brush
Battery	2 x 24V 6AH Litium ion
Charging time	5-6 hrs
Driving Distance	25km
Actual Load Capacity	120kg
Recliner	Manual 90° ~ 160°
Front / Rear Tyres	12" / 8"
Suspension Front	Perforated Front Tyre



Dinkum R120 Reclining Electric Wheelchair

Suspension Rear	2x Springs Plus Perforated Rear Tyre
-----------------	--------------------------------------

Armrest Width	450mm
---------------	-------

Cushion Width / Depth	400mm / 400mm
-----------------------	---------------

Seat to Ground	485~515mm
----------------	-----------

Widest Point	580mm
--------------	-------

Top Seat to Top Armrest	240mm
-------------------------	-------

Front to Back Length	780mm without Legrest
----------------------	-----------------------

Ground Clearance	210mm
Ground to backrest top	945mm
Unfolded Size	H595mm W590mm L500mm
Folded Size	D460mm W590mm H830mm
Joystick	360° Intelligent Brushless Display
Weight with Battery	27kg
Weight without Battery	25kg
Climbing Slope	12°
Shopper Pouch	Yes
Terrain	Semi-mid



Dinkum R120 Reclining Electric Wheelchair

Tool Kit

Yes
