



## **Etac Cloo Toilet Raiser with Armrests**



**Opening: depth x width** *29 cm x 22 cm*

---

**Seat height** *6/10/14 cm*

---

**Weight** *3.3 kg*

---

**Maximum user weight** *150 kg*

---

**Material** *Seat: Polypropylene  
Attachment:  
Aluminium  
Heightening spacers:  
Polyamide*

---

**Width between arm supports** *49 cm*

---

**Arm support height from top of seat** *21 cm*

---

**Seat width** *40 cm*

---

**Weight with arm supports** *4.8 kg*

---