



# Etac PediTurn

Tools for seated transfers between bed, wheelchair, toilet, etc.

 **etac**<sup>®</sup>

Creating Possibilities

# Etac PediTurn

## Swivel cushions for the floor

Etac PediTurn is a swivel cushion that facilitates transfer between sitting positions, such as bed to wheelchair. It is used when the user cannot move their feet during the transfer.



### Intended Use

- Support moving feet during seated transfer
- Transfer between bed and wheelchair
- Transfer between wheelchair, chair, shower chair, toilet, commode
- Transfer in and out of a car



Description	Diameter	Order Code
PediTurn	400 mm	IM99836GT



### Specifications

Top: Toughtek  
Inside: Nylon  
Underside: Toughtek  
Middle: Polyetylen

### Washing Instructions

Can be wiped off with a damp cloth.




max 150 kg  
max 330 lbs



 **etac**<sup>®</sup>  
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item