



IC333 Ultra-Lo Hospital Bed



Max SWL

Long Single: 200kg / 441 lbs SWL

King Single: 200kg / 441 lbs SWL

Double: 300kg / 661 lbs SWL

Queen: 300kg / 661 lbs SWL

Top of Bed Max Height

66 cm (25.9 inches)

Top of Bed Min Height

22 cm (8.6 inches)

Warranty

Steel Frame of Bed – 10 Year Limited Warranty

Actuators – 5 Year Limited Warranty

Hand Piece – 2 Year Limited Warranty

Control Box – 5 Year Limited Warranty

Upholstered Surround and Fabric – 3 Year Limited Warranty

Castor wheels – 5 Year Limited Warranty
