

Kmini All-Terrain Electric Wheelchair



Total length: 1150mm

Total width: 580 mm

Total height(headrest not included): 1000mm

Total weight: 55Kg

Seat to floor height(cushion not included): 440mm

Seat width: 490mm

Seat depth: 470mm

Armrest to armrest 450mm

Weight capacity:	100 kg
Turning radius:	≤0.74m
Max. driving range:	≥20 km
Max. speed:	6 km/h
Motors 250m Brushless	x 2
Footrest Adjustable	Yes
Front Caster PU Solid	8 inch
Driving wheel PU Solid	12 inch
Battery:	15.6 Ah
Ground Clearance	60mm

Height of surmountable obstacle	80mm
Width of trench crossing:	180mm
Max. slope ability:	15°
Lights at Front: LED	Yes

Easy Foldability

The Ausnew Kmini power wheelchair has the functionality to fold the wheelchair using the intelligent controller. Ideal for people that may struggle or the elderly.

Remote Control Functionality

Using a remote control device you can remotely drive your kmini into the back of your vehicle and can be quite handy when at home to steer the power chair away if it is not near you.

The headrest is fully adjustable and the seat is well contoured with a nice luxurious leather finish.

The footrest is manually adjustable to accommodate different heights or longer legs to maximise a comfortable drive for long periods.

It is worth mentioning that the kmini can be front or rear wheel drive to offer you the flexibility of driving over many terrains, obstacles and especially on roads and paves with gaps.

All in all, the Kmini folding power wheelchair is very versatile that offers you all the luxurious touches, look and feel with an affordable price tag, and is ideal for NDIS recipients and my aged care packages.

[Warranty](#)

Standard abbicare warranty – 2 years on Frame 18 months on Battery