

Leg Carriage Cushion - Relieving Leg Pillow Support



| | |
|--|---|
| TGA Listed | Yes |
| Australian Register of Therapeutic Goods Number | No. 1744454 Device: Cushion, Leg |
| Dimensions Outer | 53cm Long x 28cm Wide x 15cm High |
| Dimensions Inner | 35cm Long x 28cm Wide x 9cm High |
| Inner | 100% open cell 'breathing' urethane foam |
| Cover | Steri-Plus Outer Cover: Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester |