



Lifestyle Deluxe Self-Propelled Wheelchair



Seat Width	44 cm
Seat Depth	41 cm
Seat Height	48 cm (without cushion)
Backrest Height (from ground)	46 cm (86 cm)
Overall Height	89 cm
Overall Width	62 cm
Overall Length (with footrests)	103 cm
Armrest Height (from ground)	71 cm
Folded Height	89 cm
Folded Width	30 cm
Folded Length (with drive wheels)	62 cm (79 cm)
Total Weight (with accessories)	14.7 kg
Rear Wheels	24 inch / 60 cm Pneumatic tire
Front Wheels	15 cm Puncture proof solid tire
Material	Aluminum
SWL	140 kg
Included	Seat cushion Seatbelt Calf rest



Lifestyle Deluxe Self-Propelled Wheelchair

Backrest pad