

MyRyde Self-propelled Wheelchair



Seat Width

44 cm

Seat Depth	41 cm
Seat Height	49 cm
Backrest Height	41 cm
Overall Height	91 cm
Overall Width	62 cm
Overall Length (with footrests)	102 cm
Armrest Height (from ground)	42 cm
Folded Height	73 cm
Folded Width	33 cm
Folded Length (without footrests)	71 cm
Total Weight (with accessories)	14 kg
Rear Wheels	22 inch / 55 cm Puncture proof solid tire Lean profile for high maneuverability
Front Castors	7 inch / 17 cm Puncture proof solid tire
Material	Aluminium
SWL	120 kg
Included	Backrest pad Seat pad Seatbelt Calf rest Anti-tipping safety wheels (removable)



MyRyde Self-propelled Wheelchair