

Power SLA Electric Wheelchair - Side Folding



Seat Width

46 cm

Seat Depth

42 cm

Seat Height	52 cm
Backrest Height	41 cm
Overall Height	96 cm
Overall Width	67 cm
Overall Length (with footrests)	107 cm
Armrest Height (from ground)	73 cm
Folded Width	40 cm
Folded Length (without footrests)	70 cm
Total Weight (inc accessories & battery)	39 kg
Rear Wheels	12 inch / 30 cm Pneumatic tire
Front Castors	8 inch / 20 cm Puncture proof solid tire
Material	Steel
Range	>10 km
Max. Speed	6km/h
Motor	Two (×2) 250w with dual stage variable speed transmission