

## **Pride 130XL Pathrider Mobility Scooter**



<b>Range</b>	Up to 24 km
--------------	-------------

---

<b>Speed</b>	Up to 12 km/h
--------------	---------------

---

<b>Max SWL</b>	158 kg
----------------	--------

---

<b>Length</b>	1352 mm
---------------	---------

---

<b>Width</b>	640 mm
--------------	--------

---

<b>Heaviest Part</b>	43 kg
----------------------	-------

---

<b>brand</b>	Pride
--------------	-------

---

<b>manufacturers-part-number</b>	PR130XL-4R
----------------------------------	------------

---

<b>Tyre Type</b>	Pneumatic
------------------	-----------

---

<b>Battery Type</b>	2 X 55AH Gel
---------------------	--------------

---

<b>Motor Size (Wattage)</b>	3.09 HP – 2304W
-----------------------------	-----------------

---