

## Pride 140XL Pathrider Mobility Scooter



**Range** Up to 35 km

---

**Maximum Speed** 15km/h

---

**Backrest Height w/o Headrest** 470 mm

---

**Max SWL** 181 kg

---

**Seat Depth** 445 mm

---

**Seat Width** 495 mm

---

**Heaviest Part** 93.21 kg

---

**Brand** Pride

---

**manufacturers-part-number** PR140-R

---

**Battery Type** 2 x 75AH Gel Batteries installed

---

**Motor Size (Wattage)** 4.63 HP – 3456W

---

**Wheel Size** 13"

---

**Tyre Type** Low profile Pneumatic

---

**Maximum Incline Angle** 10 Degrees

---