



# Etac Swivel Cushions

Tools for positioning when sitting.

 **etac**<sup>®</sup>

Creating Possibilities

# AutoTurn soft swivel cushion

The Etac AutoTurn facilitates getting in and out of a car. The low friction allows the user to slide and turn into the desired position in the car seat. This minimizes strain and twisting of feet, knees, hip and back. Auto-Turn is fixed around the car seat with a belt and buckle. Handles at the sides can be used to assist the user into place, and small straps at the side can be attached around the seat belt, firstly as a location aid and then secondly as additional security.



## Intended Use

- In and out of a car
- In and out of bed



**Description**  
AutoTurn

**Diameter**  
Ø45 cm 17.7"

**Item no.**  
IM99999



## Specifications

Top: Velour (polyester)

Inside: Nylon

Underside: Nylon with PU-coating

## Washing Instructions

Can be wiped with a damp cloth

max 150 kg  
max 330 lbs



# EasyTurn soft swivel cushion

Etac EasyTurn is a soft swivel cushion that makes it easier for the user to get in and out of a car or bed. EasyTurn is portable and easy to use. It can be used in workplaces where the user often gets in and out of the car.



## Intended Use

- In and out of a car
- In and out of bed



Description	Diameter	Item no.
EasyTurn	Ø40 cm 15.7"	IM999405
EasyTurn	Ø45 cm 17.7"	IM990454
EasyTurn	Ø50 cm 19.7"	IM999505

## Specifications

Top: Velour (polyester)  
Inside: Nylon  
Backside: Nylon with polyurethan coating  
Middle: Polyetylen

## Washing Instructions

Can be wiped with a damp cloth


max 150 kg  
max 330 lbs





 **etac**<sup>®</sup>  
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item