

Bariatric Cushion



Large Size	50cm (w) x 45cm (d) x 10cm (h)
Extra Large Size	56cm (w) x 45cm (d) x 10cm (h)
Outer Cover - Dura-Fab	Woven polyester - Suede Feel
Outer Cover	Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester backing
Inner	High Density - 100% open cell 'breathing' urethane foam and visco-elastic Memory Foam.
Inner Foam Construction	2.5cm memory foam + 5cm high density max support foam + 2.5cm memory foam
Country of Manufacture	Australian Made