



## Active Ultra Lite Wheelchair



## SPECIFICATIONS

<b>Frame Material</b>	<b>Aluminum</b>	<b>Frame Material</b>
<b>Net Weight</b>	<b>8.9 kg</b>	<b>Net Weight</b>
<b>Seat Height (seat to floor)</b>	<b>50 cm</b>	<b>Seat Height (seat to floor)</b>
<b>Overall Height (floor to handles)</b>	<b>93 cm</b>	<b>Overall Height (floor to handles)</b>
<b>Seat Width</b>	<b>46 cm (18")</b>	<b>Seat Width</b>
<b>Seat Depth</b>	<b>41 cm</b>	<b>Seat Depth</b>
<b>Backrest Height (from seat)</b>	<b>44 cm</b>	<b>Backrest Height (from seat)</b>
<b>Armrest Height (from seat)</b>	<b>23-33 cm</b>	<b>Armrest Height (from seat)</b>
<b>Front Castor Size</b>	<b>8"</b>	<b>Front Castor Size</b>
<b>Rear Wheel Size</b>	<b>12"</b>	<b>Rear Wheel Size</b>
<b>Safe Work Load (SWL)</b>	<b>135 kg</b>	<b>Safe Work Load (SWL)</b>
<b>Frame Material</b>	<b>Aluminum</b>	<b>Frame Material</b>