

Auto Powered Bath Lift



SPECIFICATIONS

Overall length(with backrest in a reclined position)	985 mm
Overall width(seat incl. side flaps)	725 mm
Width of seat	400 mm
Depth of seat	500 mm
Overall length of seat	590 mm
Sitting height min	80 mm
Sitting height max	487 mm
Height of backrest	650 mm
Width of backrest	350-360 mm
Backrest adjustment	5° to max 45°
Width of baseplate	310 mm
Length of baseplate	612 mm
Load capacity max	140 kg
Total weight (excl. controller)	9.6k g
Total weight (inc. controller)	10.1 kg
Duration of lifting/lowering movement under 105kg load	Lifting:90s Lowering:50s



Auto Powered Bath Lift

Auto power off	30 min
Battery	Lithium-ion, 14.8V, 2.2 Ah
