



Comfy Lite Self-Propelled Wheelchair



SPECIFICATIONS

Type	Self-propelled	Self-propelled	Self-propelled
Frame Material	Aluminium	Aluminium	Aluminium
Net Weight	11.7 kg	11.9 kg	12.1 kg
Seat Height (seat to floor)	45 cm	45 cm	45 cm
Overall Height (floor to handles)	92.5 cm	92.5 cm	92.5 cm
Seat Width	40 cm	45 cm	47.5 cm
Seat Depth	42 cm	42 cm	42 cm
Backrest Height (from seat)	42.5 cm	42.5 cm	42.5 cm
Armrest Height (from seat)	23 cm	23 cm	23 cm
Front Castor Size	6"	6"	6"
Rear Wheel Size	22"	22"	22"
Safe Work Load (SWL)	135 kg	135 kg	135 kg