

Euro Lite Seat Walker



SPECIFICATIONS

Frame Material	Aluminium	Carbon Fibre	Aluminium	Carbon Fibre	Aluminium	Carbon Fibre
Size (height)	Low	Low	Standard	Standard	Tall	Tall
Product Frame Weight	6.4 kg	5.2 kg	6.6 kg	5.3 kg	6.8 kg	5.4 kg
Product Gross Weight	7.8 kg	6.6 kg	8 kg	6.7 kg	8.2 kg	6.8 kg
Seat Height (seat to floor)	50 cm	50 cm	57 cm	57 cm	63 cm	63 cm
Overall Height (floor to handles)	66 cm – 84 cm	66 cm – 84 cm	73 cm – 90 cm	73 cm – 90 cm	77 cm – 96 cm	77 cm – 96 cm
Seat Width (in between arms)	45 cm	45 cm	45 cm	45 cm	45 cm	45 cm
Outer Width (including rear castors)	59 cm	60.5 cm	59 cm	60.5 cm	59 cm	60.5 cm
Front Castor Size	9"	9"	9"	9"	9"	9"
Rear Castor Size	9"	9"	9"	9"	9"	9"
Safe Work Load (SWL)	130 kg	130 kg	130 kg	130 kg	130 kg	130 kg