

Lifestyle Extra - Self-Propelled Wheelchair



Seat Width

50 cm

Seat Depth

42 cm

Seat Height

49 cm (without cushion)

Backrest Height (from ground)	44 cm (86 cm)
Overall Height	93 cm
Overall Width	67 cm
Overall Length (with footrests)	105 cm
Armrest Height (from ground)	69 cm
Folded Height	80 cm
Folded Width	29 cm
Folded Length	81 cm
Total Weight (with accessories)	15 kg
Rear Wheels	24 inch / 60 cm Puncture proof solid tire
Front Castors	15 cm Puncture proof solid tire
Material	Aluminium
SWL	140 kg
Included	Seat cushion Seatbelt Calf rest on each footplate