

MyRide Kids - Pediatric Wheelchair



Seat Width

35 cm

Seat Depth	35 cm
Seat Height	44 cm
Backrest Height	40 cm (74 cm from ground)
Overall Height	90 cm
Overall Width	54 cm
Overall Length (with footrests)	87 cm
Armrest Height (from ground)	60 – 65 cm
Folded Height	60 cm
Folded Width	32 cm
Folded Length (without footrests)	72 cm
Total Weight (with accessories)	12.5 kg
Rear Wheels	20 inch / 50 cm Puncture proof solid tire
Front Castors	4 inch / 10 cm Puncture proof solid tire
Material	Aluminium
SWL	80 kg
Included	Backrest cushion Ergonomic seat cushion Body harness Calf rest Anti-tipping safety wheels