



MyRide Power Wheelchair - Front Folding



SPECIFICATIONS

Seat Width	42 cm
Seat Depth	43 cm
Seat Height	43 cm
Backrest Height	44 cm
Overall Height	90 cm
Overall Width	58 cm
Overall Length (with footrests)	98 cm
Armrest Height (from ground)	87 cm
Folded Width	tba cm
Folded Length (without footrests)	tba cm
Total Wight (inc accessories & battery)	21 kg
Rear Wheels	12 inch / 30 cm Pneumatic tire
Front Castors	8 inch / 20 cm Puncture proof solid tire
Material	Aluminium
Range	20 km
Max. Speed	6km/h



MyRide Power Wheelchair - Front Folding

Motor	Two (x2) 130w
Battery	24V, 12Ah Lithium
SWL	100 kg
Included	Backrest pad Seat cushion Seatbelt Anti-tipping safety wheels (removable) Charger Battery pack