

## Solax Mobifree



# SPECIFICATIONS

**Adjustable Tiller Angle**

---

Yes

**Adjustable Tiller Height**

---

Yes

**Anti-Tip Mechanism**

---

## Rear Anti-Tip Wheels

---

### Armrests

---

Yes

---

### Battery

---

24V/10ah Lithium

---

### Battery Range

---

Up to 15km

---

### Battery Weight

---

1.8kg

---

### Dimensions

---

969x530x950

---

### Dimensions when folded

---

343x530

---

### Drive Wheel

---

170mm

---

### Ground Clearance

---

40mm

---

### Maximum Speed

---

6km/h

---

### Number of Components

---

1 + Battery

---

### Slope ability

---

12°

---

### Tyres

---

Solid Puncture Proof

---

### Weight Capacity

---

125kg

---



## Solax Mobifree

**Weight with Battery**

---

22.3kg

**Weight without Battery**

---

20.5kg

---