



Super Budget Wheelchair



SPECIFICATIONS

Type	Transit	Self-Propelled
Frame Material	Steel	Steel
Net Weight	15.4 kg	17.6 kg
Seat Height (seat to floor)	48 cm	48 cm
Overall Height (floor to handles)	92 cm	92 cm
Seat Width	45 cm (18")	45 cm (18")
Seat Depth	40 cm (16")	40 cm (16")
Backrest Height (from seat)	42.5 cm	42.5 cm
Armrest Height (from seat)	26 cm	26 cm
Front Castor Size	8"	8"
Rear Wheel Size	12.5"	24"
Safe Work Load (SWL)	115 kg	115 kg