

Bathroom

Shower Stool w Padded Seat

The Affinity Shower Stool with soft PU seat offers the ultimate in comfort. Strong sturdy construction gives this chair an exceptional weight rating of 150kg. The soft removable PU seat pad makes cleaning a breeze. U shape cutout seat front assists with personal hygiene. Available in 3 seat widths and optional removable backrest.

- * Lightweight powder coated aluminium frame
- * Removable soft PU seat - easy cleaning
- * Drain holes to reduce slip and water build up
- * Height adjustable—5 choices at 2.5cm increments

- * Large non-slip stoppers
- * Available in 3 widths
- * Optional extra padded backrest
- * Weight capacity 150kg












Product Code: SC42..... 420mm
 SC47..... 470mm
 SC51..... 510mm



Product Code: SC51BR

Specifications

overall height	overall length	overall width	seat depth	seating width	seat height	armrest height	unit weight	weight capacity
								
42 60-68	54	51	33	42	40-48	20	3.6	150
47 60-68	54	56	33	47	40-48	20	4.1	150
51 60-68	54	60	33	51	40-48	20	5	150

Measurements in cm, Weight in kg