

Seating

Astral Perching Stool

The Astral Perching Stool has been developed for comfort and longevity, the contoured lines of both the seat and back are not only attractive to the eye, but also offer excellent ergonomic support with seamless polyurethane seat. The height adjustable frame has side mounted stress bars and a large footprint for safety and stability. Ideal for those who struggle to stand for long periods.

- * Ideal for perching when working in the kitchen or at a bench
- * Adjustable seat height
- * Padded backrest
- * Easy clean
- * Weight Capacity 159kg



Product Code: VG863



Product Code: VG865

Specifications

Seat Width 350mm
Seat Depth 300mm
Seat Height 500-650mm
Width Between Arms 480mm
Arm Height 170mm
Back Height 270mm
Footprint 410-460 x 585-650mm
Weight Capacity 159kg