

## Bathroom *Bariatric*

### Bath Transfer Bench

The durable aluminium Bath Transfer Bench is placed over the bath for those who have a shower over the bath situation. The wider seat allows the client to slide across from the outside of the bath. Two large suction feet for secure placement. Fully height adjustable will fit over most baths.

Optional extension legs available for external leg for higher bath tubs.



Product Code: RM424

- \* Durable aluminium frame
- \* Detachable back rest, can be attached to either side
- \* Heavy duty suction tips
- \* Detachable arm rest
- \* Weight capacity 200kg

- \* Optional extension legs available to raise bench into higher position

- **Complies with Australian Standards  
AS/NZS 3973:2009**

#### Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Tub Clearance Height	Armrest Height	Unit Weight	Weight Capacity
84-99	52	78	41	68	41-55	30	6	200

Measurements in cm, Weight in kg

### Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com