

Bathroom Bariatric

Bath Transfer Bench HDuty

The heavy duty durable aluminium Bath Transfer Bench is placed over the bath for those who have a shower over the bath situation. The wider seat allows the client to slide across from the outside of the bath. Two large suction feet for secure placement. Fully height adjustable will fit over most baths.

Optional extension legs available for external leg for higher bath tubs.

- * Durable aluminium frame
- * Moulded seat and removeable backrest
- * Heavy duty suction tips
- * Detachable armrest
- * Weight capacity RM424 180kg
- * Weight capacity RM425 280kg

- * Optional extension legs available to raise bench into higher position

- **Complies with Australian Standards
AS/NZS 3973:2009**



Specifications

	overall height	overall depth	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
RM424	85-97.5	58	81	46	66	45-55	20	4.7	180
RM425	83-96	58	79	41.5	68	46-59	26	6.6	280

Measurements in cm, Weight in kg