


## Bedroom



### Bed Rails

Place under mattress to assist user in getting in and out of bed.



#### Specifications


Dimensions (mm).....520x810  
 Height (mm).....440-570  
 Weight capacity (kg).....127

- \* Height adjustable safety rail
- \* Large base for improved stability
- \* High quality easy to clean epoxy resin finish

Product Code: RM506-KD

#### Specifications

Dimensions (mm).....580x600x680  
 Height (mm).....580-700  
 Weight capacity (kg).....127



Product Code: VY428



Product Code: VY428S

- \* Designed for use on slatted beds with a lip
- \* Large extendable base for improved stability

#### Specifications

Overall dimensions (mm).....500x450x580-820  
 Handle Width ..... 330  
 Weight capacity (kg).....160

The double bedrail passes below the mattress and the looped handles extend up either side of the bed.

- \* Easy adjustable telescopic tube to adjust width
- \* One side can be turned down if not required on either sides of the bed
- \* High quality, durable chrome finish

#### Specifications

Width Adjustable (mm)....910-1530  
 Height (mm) ..... 360  
 Depth (mm) ..... 340  
 Weight Capacity (kg) .....127



Product Code: VY430