

NEW



Invacare® **Comet™ Alpine**

Powerful, Determined, Adventurous: Go Far with the Invacare Comet Alpine

New to the Invacare Comet family, the **Comet Alpine** version offers improved traction and a larger motor, ideal for challenging outdoor uses. The robust power unit and 4-point independent suspension system ensures a continuous comfortable ride even on bumpy ground and can tackle the most challenging of hills. The **Comet Alpine** includes a variety of safety features such as speed reduction when turning, brake light when stopping and a 2-step disengaging lever to prevent the scooter from freewheeling.



➤ Extra Power

Beyond its stylish appearance, the **Comet Alpine** proves itself a real winner. The scooter will win you over with its powerful driving characteristics and impressive maneuverability going uphill and even when driving downhill or on loose gravel. The stability of the **Comet Alpine** in challenging surroundings leaves no doubt that this scooter is out of the ordinary.

➤ Superior Range of Safety Features

The Automatic Speed Reduction feature automatically reduces your speed around bends, preventing you from going too fast when turning.

➤ Challenge Your Comet

With the **Comet Alpine** you can leave limitations behind, enjoying outings across uneven or wet ground. The **Comet Alpine** has built-in anti-splash guards, which protect vital components from water and dirt.



Automatic Speed Reduction

Features & Accessories



Seat suspension module*
Adds extra comfort during driving.



Seat Raiser*
Electronically raises the seat to the preferred height.



Off-board charging
An off-board charger kit is available, which allows you to conveniently charge the batteries away from the scooter.

*Only available for Comet and Comet Alpine.

Technical Data

								Model: 4 Wheel
51cm	48cm	44-51.5cm	50cm	66cm	147.5cm	143kg*	160kg	Colours: Blue, Silver & Red
								Wheels: Front & Rear 31cm/13" Pneumatic
(2) 12V 75Amp	163cm	10cm	13.5cm	10°	Cont: 650W Peak: 1800W	50km	10 km/h	Suspension: Front & Rear Spring & Adjustable
								Seat Type: Extra padding, ergonomic, recliner lever, and headrest
								Std. Batt. Charger: 24V / 8 Amp (off-board)
								Electronics: 160 Amp Rhino Controller

Compliance: TGA-registered (ARTG# 119338) CE Marked TUV certified – tested to European Standard (DIN EN 12184:1999)

*Without batteries